

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41T01-85

April 1, 2005

[3rd edition]

Product Name: TYRAMINE OXIDASE [TOD]

(Amine:oxygen oxidoreductase(deaminating), EC 1. 4. 3. 4)

Product Number: T-25

Origin: *Arthrobacter* sp

Description: White to light brownish lyophilized powder

Activity: ≥ 3 U/mg

Appearance of solution:
(5U/ml in buffer solution)


| |
|---------------------|
| Clear |
| A400nm ≤ 0.150 |
| A500nm ≤ 0.060 |
| A650nm ≤ 0.030 |

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant