



POINT-OF-CARE



OSOM® BVBLUE® Test

FOR THE DETECTION OF SIALIDASE ACTIVITY FROM GARDNERELLA, BACTEROIDES, PREVOTELLA, & MOBILUNCUS

Features and Benefits


- **Accurate** – Test and treat bacterial vaginosis with confidence
- **Fast** – Less than one minute of hands on time. Results in 10 minutes
- **Objective** – Instant color change provides clear, easy-to-read results
- **Exceptional Customer Service** – Online training modules and technical support staffed by certified medical technologists/professionals

Performance Characteristics


Sensitivity/Specificity	Sample Type	Limit of Detection	Storage Temperature	Time to Result
92.8%/98.0%*	Vaginal swab	7.64 U**	2-8°C/36-46°F	10 minutes

*Compared to gram stain


** Mean vaginal sialidase activity level in healthy controls is 2.7 U. See IFU for additional information



OSOM® BVBLUE® Test



SIMPLE TO USE / RESULTS IN MINUTES



Ordering Information

Description	Configuration	Catalog Number	Storage	CPT Code
OSOM BVBLUE TEST	25 Tests	183	Refrigerated	87905QW
OSOM BVBLUE CONTROL KIT	2 X 5 ml (Pos and Neg)	184	Refrigerated	
OSOM VAGINITIS TRANSPORT SET (TRICH & BV)	50 Swab	185	Room Temp	
OSOM BVBLUE TRANSPORT SWAB	30 Swabs	186	Room Temp	

To learn more about all of our OSOM Tests scan the code or go to osomtests.com



THE AMERICAS

SEKISUI Diagnostics, LLC
 One Wall Street
 Burlington, MA 01803
 Phone: 800 332 1042
 Fax: 800 762 6311
 Email: questions@sekisui-dx.com



Because every result matters™

sekisuidiagnostics.com

