POINT-OF-CARE



OSOM[®] Ultra Strep A Test

STRONGER CLINICAL PERFORMANCE TAKES LATERAL FLOW TESTING TO THE NEXT LEVEL

Features and Benefits

- Accurate Statistically equivalent to single swab culture*
- Dual Labeling Technology Greater sensitivity than conventional single labeling technology**
- Fast and Easy workflow Pre-filled ampules (cat # 147) -two minutes of hands on time
- Simple CLIA-waived
- Exceptional Customer Service Online training modules and technical support staffed by certified medical technologists/professionals

Performance Characteristics						
Sensitivity/Specificity	Sample Type	Limit of Detection	Storage Temperature	Time to Result		
96.4% / 96.3%***	Throat Swab	<10 colonies	15-30°C / 59-86°F	7 Minutes		

*Compared to rigourous gold standard (RGS). Refer to IFU


**Refer to IFU

***Compared to single swab culture






OSOM[®] Ultra Strep A Test



INTERPRETATION OF RESULTS



Ordering Information						
	Configuration	Catalog Number	CPT Code (US only)			
OSOM ULTRA STREP A TEST KIT/PRE-FILLED VIALS	25 tests	147	87880QW			
OSOM ULTRA STREP A TEST KIT/DROPPER VIALS	50 tests	149				

To learn more about all of our OSOM Tests scan the code or go to <u>osomtests.com</u>

SEKISUI Diagnostics, LLC One Wall Street Burlington, MA 01803 Phone: 800 332 1042 Fax: 800 762 6311 Email: questions@sekisui-dx.com





© 2021 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUi Diagnostics, LLC.

