



SPECIFICATION

ISO: 75F41P16-85
June 9, 2017
【 1st edition 】

Product Name:	PHOSPHOLIPASE-D [PLD II] (Phosphatidylcholine phosphatidohydrolase, EC 3. 1. 4. 4)
Product Number:	T-222
Origin:	<i>Streptomyces chromofuscus</i>
Description:	Light purple to brownish lyophilized powder
Activity:	30~90 U/mg
Appearance of solution: (10U/mL in buffer solution)	Clear A400nm ≤0.200 A500nm ≤0.060 A650nm ≤0.030
Contaminating enzymes:	
Catalase	≤0.6 % (U/U)
Glucose oxidase	≤0.02 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

S.ISHIKAWA
Plant Director
Ohito diagnostics plant

E.OHKI
Quality assurance and Quality control Manager
Ohito diagnostics plant