



SPECIFICATION

ISO: 75F41H06-85
September 17, 2010
【 1st edition 】

Product Name: 12 α -HYDROXYSTEROID DEHYDROGENASE [12 α -HSD II]
(12 α -Hydroxysteroid:NAD 12 α -oxidoreductase, EC 1. 1. 1. 176)

Product Number: T-190
Origin: *Bacillus sphaericus*
Description: White to light brownish lyophilized powder
Activity: ≥ 100 U/mg

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant