

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41H02-85

April 1, 2005

[3rd edition]


Product Name: 3-HYDROXYBUTYRATE DEHYDROGENASE [3-HBDH II]

(D-3-Hydroxybutyrate:NAD oxidoreductase, EC 1.1.1.30)


Product Number:	T-68
Origin:	<i>Alcaligenes faecalis</i>
Description:	White lyophilized powder
Activity:	≥ 1500 U/mg
Appearance of solution: (1mg/ml in buffer solution)	Clear A400nm ≤ 0.050 A500nm ≤ 0.040 A650nm ≤ 0.030
Contaminating enzymes: NADHoxidase	$\leq 0.0004\%$ (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant