

Asahi Kasei Pharma Corporation  
Fine Chemicals & Diagnostics Division  
9-1 Kanda Mitoshiro-cho  
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department  
Phone +81-(0)3-3259-5775  
Fax +81-(0)3-3259-5741  
diagnostics@om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41A03-85

April 1, 2005

[ 3<sup>rd</sup> edition ]

Product Name: ACYL-COA SYNTHETASE [ACS]

(Acid:CoA Ligase(AMP), EC6. 2. 1. 3)

Product Number: T-16

Origin: *Pseudomonas fragi*

Description: White lyophilized powder

Activity: 2~8 U/mg

Appearance of solution: Clear  
(2U/ml in buffer solution )  
A400nm  $\leq 0.050$   
A500nm  $\leq 0.025$   
A650nm  $\leq 0.012$


Contaminating enzymes:  
Catalase  $\leq 1.0$  (U/U)  
Glucose oxidase  $\leq 0.1\%$  (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*

*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T.OTANI  
Plant Director  
Ohito diagnostics plant



H.NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant