

Asahi Kasei Pharma Corporation  
Fine Chemicals & Diagnostics Division  
9-1 Kanda Mitoshiro-cho  
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department  
Phone +81-(0)3-3259-5775  
Fax +81-(0)3-3259-5741  
diagnostics@om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41A09-85

April 1, 2005

[ 3<sup>rd</sup> edition ]

Product Name: ALKALINE PHOSPHATASE [ALP]

(Orthophosphoric-monoester phosphohydrolase , EC 3. 1. 3. 1)

Product Number: T-08

Origin: *Escherichia coli*

Description: White lyophilized powder

Activity:  $\geq 40$  U/mg

Appearance of solution: Clear  
(50U/ml in buffer solution)


A400nm	$\leq 0.015$
A500nm	$\leq 0.015$
A650nm	$\leq 0.015$

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*

*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T. OTANI  
Plant Director  
Ohito diagnostics plant



H. NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant