

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41A10-85
April 1, 2005
[3rd edition]

Product Name: ASCORBATE OXIDASE [ASOM]
(L-Ascorbate:oxygen oxidoreductase , EC 1. 10. 3. 3)

Product Number: T-53

Origin: *Acremonium* sp.

Description: Light blue lyophilized powder

Activity: 200~450 U/mg


Appearance of solution: Clear
(50U/ml in buffer solution) A400nm \leq 0.010
A500nm \leq 0.010
A650nm \leq 0.010

Contaminating enzymes:
Catalase \leq 0.02% (U/U)
ATPase \leq 0.001% (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant