

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41C14-85

April 1, 2005

[3rd edition]

Product Name: CHOLESTEROL OXIDASE [CON II -FD]

(Cholesterol : oxygen oxidoreductase, EC 1. 1. 3. 6)

Product Number: T-101

Origin: *Microorganism*

Description: Yellowish lyophilized powder

Activity: ≥ 15 U/mg

Appearance of solution: Clear
(5U/ml in buffer solution)
A400nm ≤ 0.100
A500nm ≤ 0.050
A650nm ≤ 0.010


Contaminating enzymes:
Catalase ≤ 0.5 % (U/U)
Cholesterol esterase ≤ 0.02 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant