

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41L09-85

April 1, 2005
[3rd edition]

Product Name: LIPASE [LPM]
(Triacyl glycerol acylhydrolase, EC 3. 1. 1. 3)

Product Number: T-116

Origin: *Chromobacterium viscosum*

Description: White to pale brownish lyophilized powder

Activity: 1500~3500 U/mg

Appearance of solution: Clear
(1mg/ml in buffer solution)
A400nm \leq 0.100
A500nm \leq 0.050
A650nm \leq 0.050

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant