

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41L10-85
November 18, 2005
【 2nd edition 】

Product Name: LYSOPOSPHOLIPASE [LYPL]

 (2-Lysophosphatidylcholine acylhydrolase, EC 3. 1. 1.5)

Product Number: T-32

Origin: *Vibrio* sp.

Description: White to brownish lyophilized powder

Activity: ≥ 20.0 U/mg

Appearance of solution: Clear
(1mg/ml in buffer solution) A400nm ≤ 0.200
 A500nm ≤ 0.100
 A650nm ≤ 0.050

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.