

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41P08-85
April 1, 2005
[3rd edition]


Product Name: PURINE-NUCLEOSIDE PHOSPHORYLASE [PNPL II]

(Purine-nucleoside phosphorylase, EC 2.4.2.1)

Product Number:	T-69
Origin:	<i>Bacillus</i> sp.
Description:	Colorless to light brownish solution
Activity:	≥ 500 U/ml
Appearance of solution: (10U/ml in buffer solution)	Clear A400nm ≤ 0.010 A500nm ≤ 0.010 A650nm ≤ 0.010
Contaminating enzymes: NADHoxidase	$\leq 0.002\%$ (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant