

# Glucose Dehydrogenase PQQ

Catalogue No. GD-2  
E.C. Number 1.1.99.17

## ORIGIN

Bacterial

## KEY FEATURES

- o RECOMBINANT
- o HIGH ACTIVITY POWDER (> 8000 U/MG)
- o APOFORM, REQUIRES PQQ FOR ACTIVATION

## ASSAY PRINCIPLE

Reacts with glucose (with suitable acceptor) to form glucono 1,5 lactone.

## APPLICATION

Can be used as a signal amplifier.

## PRODUCT FORMAT

Freeze-dried powder

**DP** DISCOVERY PRODUCTS ARE ENZYMES FROM OUR R&D LABS PRIOR TO FULL COMMERCIALISATION



**EMEA**  
Sekisui Diagnostics (UK) Ltd  
50 Gibson Drive,  
Kings Hill, West Malling,  
Kent, ME19 4AF, UK  
Phone: +44 (0)1732 220022  
Fax: + 44 (0)1732 220024

**AMERICAS**  
Sekisui Diagnostics, LLC  
4 Hartwell Place  
Lexington,  
MA 02421, USA  
Phone: 800 332 1042  
Fax: 800 762 6311

**JAPAN AND ASIA GROUP**  
Sekisui Medical Co., Ltd.  
13-5, Nihonbashi 3-Chome,  
Chuo-ku, Tokyo 103-0027 Japan  
Phone: +81-3-3272-0828  
Fax: +81-3-3272-0907

Email: [engage@sekisuienzymes.com](mailto:engage@sekisuienzymes.com)