

Technopath Multichem® S Plus Control Values for Acetaminophen

ACETAMINOPHEN VALUES						
Multichem S Plus Control	N	Mean (µg/mL)	Standard Deviation (SD)	3SD	20% Mean	Established Range (Mean ±20%)
Level 1 Lot 18004181, Exp. 2020-10-31	32	14.85	0.18	0.54	2.97	11.88 – 17.82
Level 2 Lot 18005182, Exp. 2020-10-31	32	68.88	0.41	1.23	13.78	55.10 – 82.66
Level 3 Lot 18005183, Exp. 2020-11-30	32	211.92	1.20	3.60	42.38	169.54 – 254.30

Values are provided by Sekisui Diagnostics for use with Sekisui Acetaminophen 506 reagent (Abbott List Number 03R11-20) as determined using the Abbott Architect c Systems. Please follow all instructions for use by the analyzer and controls manufacturers.

Effective Date: 11-Mar-2019

1. Multichem is a registered trademark of Technopath Clinical Diagnostics, Tipperary, Ireland.
2. Architect is a trademark of Abbott Laboratories in various jurisdictions.