

## Technopath Multichem® S Plus Control Values for Acetaminophen

<b>Multichem S Plus Control</b>	<b>N</b>	<b>Mean (µg/mL)</b>	<b>Standard Deviation (SD)</b>	<b>3SD</b>	<b>20% Mean</b>	<b>Established Range (Mean ±20%)</b>
Level 1 Lot 18903191, Exp. 2021-08-31	32	7.82	0.15	0.45	1.56	<b>6.26 – 9.38</b>
Level 2 Lot 18903192, Exp. 2021-08-31	32	64.65	0.41	1.23	12.93	<b>51.72 – 77.58</b>
Level 3 Lot 18903193, Exp. 2021-08-31	32	198.09	0.77	2.31	39.62	<b>158.47 – 237.71</b>

Values are provided by Sekisui Diagnostics for use with Sekisui Acetaminophen 506 reagent (Abbott List Number 03R11-20) as determined using the Abbott Architect c Systems. Please follow all instructions for use by the analyzer and controls manufacturers.

Effective Date: 07-Feb-2020

1. Multichem is a registered trademark of Technopath Clinical Diagnostics, Tipperary, Ireland.
2. Architect is a trademark of Abbott Laboratories in various jurisdictions.