

Technopath Multichem® S Plus Control Values for Acetaminophen

ACETAMINOPHEN VALUES						
Multichem S Plus Control	N	Mean (µg/mL)	Standard Deviation (SD)	3SD	20% Mean	Established Range (Mean ±20%)
Level 1 Lot 17610171, Exp. 2020-03-31	32	14.56	0.19	0.57	2.91	11.65 – 17.47
Level 2 Lot 17610172, Exp. 2020-03-31	32	71.32	1.43	4.29	14.26	57.06 – 85.58
Level 3 Lot 17610173, Exp. 2020-03-31	32	221.95	6.32	18.96	44.39	177.56 – 266.34

Values are provided by Sekisui Diagnostics for use with Sekisui Acetaminophen 506 reagent (Abbott List Number 03R11-20) as determined using the Abbott Architect c Systems. Please follow all instructions for use by the analyzer and controls manufacturers.

Effective Date: 08-Jun-2018

1. Multichem is a registered trademark of Technopath Clinical Diagnostics, Tipperary, Ireland.
2. Architect is a trademark of Abbott Laboratories in various jurisdictions.