

Technopath Multichem® S Plus Control Values for Acetaminophen

| ACETAMINOPHEN VALUES | | | | | | |
|--|----|--------------|-------------------------|------|----------|-------------------------------|
| Multichem S Plus Control | N | Mean (µg/mL) | Standard Deviation (SD) | 3SD | 20% Mean | Established Range (Mean ±20%) |
| Level 1 Lot 17205171, Exp. 2019-11-30 | 32 | 14.00 | 0.36 | 1.08 | 2.80 | 11.20 – 16.80 |
| Level 2 Lot 17205172, Exp. 2019-11-30 | 32 | 67.22 | 1.13 | 3.39 | 13.44 | 53.78 – 80.66 |
| Level 3 Lot 17205173, Exp. 2019-11-30 | 32 | 203.36 | 3.00 | 9.00 | 40.67 | 162.69 – 244.03 |

Values are provided by Sekisui Diagnostics for use with Sekisui Acetaminophen 506 reagent (Abbott List Number 03R11-20) as determined using the Abbott Architect c Systems. Please follow all instructions for use by the analyzer and controls manufacturers.

Effective Date: 22-Jan-2018

1. Multichem is a registered trademark of Technopath Clinical Diagnostics, Tipperary, Ireland.
2. Architect is a trademark of Abbott Laboratories in various jurisdictions.