

What is QC Inside®?



HOW YOU CAN IMPROVE COST EFFICIENCY WHEN PERFORMING POINT-OF-CARE TESTING

SEKISUI Diagnostics offers QC Inside® for most of our point-of-care testing products, meaning a kit of 25 tests includes 2 additional tests that can be used for external quality control testing.

At no additional charge, each kit contains:

- extra test strips
- extra extraction buffers
- or other components

... allowing for an additional 2 tests per kit.

Learn more about QC Inside® by visiting sekisuidiagnostics.com/qc-inside/



Why we developed QC Inside®?

For all CLIA-waived tests, it is recommended to test one positive and one negative control every time there is a:

- new product shipment
- change in lot number
- new untrained user

Without QC Inside®, you would have:

25 Tests per Kit


The kit configuration for most rapid point-of-care test kits.

-2 External Quality Control Tests

The number of tests consumed when QC testing is deemed necessary.

=23 Patient Tests per Kit

The number of patient tests leftover per kit when QC testing is concluded.





Which SEKISUI Diagnostics point-of-care tests have QC Inside®?

Product Name	Catalog Number	Tests per Kit	Contains QC Inside®?
Respiratory Health			
OSOM® Flu SARS-CoV-2 Combo Test	1080	25	Yes
Acucy® Influenza A&B Test	1010	25	Yes
OSOM® Ultra Plus Flu A&B Test	1032	25	Yes
OSOM® Mono Test	145	25	Yes
OSOM® COVID-19 Antigen Rapid Test	1066-40	40	No
OSOM® Strep A Test	141	50	Yes
OSOM® Ultra Strep A Test	149	50	Yes
OSOM® Ultra Strep A Test	147	25	No
Women's Health			
OSOM® BVBLUE® Test	183	25	Yes
OSOM® Trichomonas Test	181	25	No
OSOM® Ultra hCG Combo Test	1004	25	Yes
OSOM® hCG Card Pregnancy Test	102	25	No
Gastrointestinal Health			
OSOM® iFOB Test	1002	25	Yes
OSOM® H. Pylori Test	175	25	No



© 2024 SEKISUI Diagnostics, LLC. All rights reserved. Acucy®, OSOM® and QC Inside® are registered trademarks of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUI Diagnostics, LLC.

SEKISUI
DIAGNOSTICS

Because every result matters™

80-8952-00-01 04/24